

Phaco Tips

You can use the following products with the **WHITESTAR Signature** System, the **SOVEREIGN** and the **SOVEREIGN Compact** System.

Disposable **LAMINAR** Flow Phaco Tips:

21 Gauge

Part #	Name	Quantity
OPOR1521G	Tip, 15-degree, with infusion sleeve	6 per box
OPOR3021G	Tip, 30-degree, with infusion sleeve	6 per box
OPOCR3021L	Tip, 30-degree, curved, with infusion sleeve	6 per box

20 Gauge

Part #	Name	Quantity
OPOCR3020L	Tip, 30-degree, curved, with infusion sleeve	6 per box
OPOD3020L	Tip, 30-degree, with infusion sleeve	6 per box

19 Gauge

Part #	Name	Quantity
OPOD3019L	Tip, 30-degree, with infusion sleeve	6 per box
OPOF3019L	Tip, 30-degree, flared, with infusion sleeve	6 per box

For use with **OPOMTWL** metal tip wrench



Metal Tip Wrench
Part # **OPOMTWL**

Multi-Use **LAMINAR** Flow Phaco Tips:

20 Gauge

Part #	Name	Quantity
OPOR1520L	Tip, 15-degree, with infusion sleeve	6 per box
OPOR3020L	Tip, 30-degree, with infusion sleeve	6 per box

19 Gauge

Part #	Name	Quantity
OPOR0019L	Tip, 0-degree, with infusion sleeve	6 per box
OPOR1519L	Tip, 15-degree, with infusion sleeve	6 per box
OPOR3019L	Tip, 30-degree, with infusion sleeve	6 per box

Sleeves, test chambers and miscellaneous:

Part #	Name	Quantity
OPOS19L	Multi-use, 19-gauge, LAMINAR sleeve and test chamber	6 per box
OPOS20L	Multi-use, 20-gauge, LAMINAR sleeve and test chamber	6 per box
OPOS21LD	Disposable, 21-gauge, LAMINAR sleeve and test chamber	6 per box
OPOMTWL	Metal tip wrench	1 each



Disposable **LAMINAR Sleeve and Test Chamber, 21-gauge**
Part # **OPOS21LD**






Multi-use **LAMINAR Sleeve and Test Chamber, 20-gauge**
Part # **OPOS20L**



Multi-use **LAMINAR Sleeve and Test Chamber, 19-gauge**
Part # **OPOS19L**

Phaco Tips

Options of AMO Phaco tips & sleeves combinations

	Incision size 2.2 - 2.4 mm	Incision size 2.4 - 2.7 mm	Incision size 2.7 - 3.0 mm
Sleeve	light blue (OPOS21LD) 	yellow (OPOS20L) 	dark blue (OPOS19L) 
Tips	21g tips OPOR1521G, disposable OPOR3021G, disposable OPOCR3021L, disposable	21g tips OPOR1521G, disposable OPOR3021G, disposable OPOCR3021L, disposable	—
	—	20g tips OPOR1520L, reusable OPOR3020L, reusable OPOD3020L, disposable OPOCR3020L, disposable	—
	—	—	19g tips OPOR0019L, reusable OPOR1519L, reusable OPOR3019L, reusable OPOD3019L, disposable OPOF3019L, disposable



0° Bevel



15° Bevel



30° Bevel



45° Bevel



Straight Tip



Bent Tip
20° Angle



LAMINAR Flow Phaco Tip – Straight



LAMINAR Flow Phaco Tip – Curved



LAMINAR Flow Phaco Tip – Flared

Irrigation and Aspiration Accessories

You can use the following products with all **AMO** Phaco Systems.

Coaxial:

PHACOFIT Irrigation / Aspiration Products

Part #	Name	Quantity
OM055002	Titanium handle	1 each
OM0551012	Tip, curved, 0.3 mm port	1 each
OM0551013	Tip, 45 degree, 0.3 mm port	1 each
OM0551014	Tip, 90-degree, 0.3 mm port	1 each
OM0551016	Tip, straight, 0.3 mm port for silicone sleeve	1 each
OM05510116	Tip, 45-degree, 0.3 mm port for silicone sleeve	1 each
Ti-98431	Tip, 21-gauge, straight, 0.3 mm port	1 each
Ti-98432	Tip, 21-gauge, bent, 0.3 mm port	1 each
OM05510112	Replacement O-ring	5 per box

NEW



PHACOFIT Multi-tip Irrigation and Aspiration Products

Bimanual:

Duet Bimanual 22g Irrigation / Aspiration System



NEW

Part #	Name	Quantity
DU02100	Duet aspiration handle	1 each
DU02200	Duet irrigation handle	1 each
DU-02355	Aspiration tip, 22-gauge, curved, closed end, 0.3 mm port, capsule polisher	1 each
DU-02354	Irrigation tip, 22-gauge, curved, closed end	1 each



Irrigation and Aspiration Accessories

Use of AMO irrigation & aspiration systems according to incision size

		Incision size 2.2 - 2.4 mm	Incision size 2.4 - 2.7 mm	Incision size 2.7 - 3.0 mm
I/A handles & tips	Duet Bi-Manual I/A System	MST Duet aspiration handle: DU02100 WITH	–	–
		Curved aspiration tip, closed end, 22g (0.7 mm), 0.3 mm port capsule polish: DU-02355	–	–
		AND		
		MST Duet irrigation handle: DU02200 WITH	–	–
		Curved irrigation tip, closed end, 22g (0.7 mm): DU-02354	–	–
	PhacoFit Coaxial I/A System	I/A Titanium handle: OM055002 WITH	I/A Titanium handle: OM055002 WITH	I/A Titanium handle: OM055002 WITH
		MST Barrett V2, 21g, straight, capsule polish, 0.3 mm port: Ti-98431 OR	Tip, curved, 0.3 mm port: OM0551012 OR	Tip, curved, 0.3 mm port: OM0551012 OR
		MST Barrett V2, 21g, bent, capsule polish, 0.3 mm port: Ti-98432	Tip, 45-degree, 0.3 mm port: OM0551013 OR	Tip, 45-degree, 0.3 mm port: OM0551013 OR
		OR:	Tip, 90-degree, 0.3 mm port: OM0551014 OR	Tip, 90-degree, 0.3 mm port: OM0551014 OR
		Tip, 45-degree/Tip, straight: OM05510116/OM0551016 WITH 21g sleeve OPOS21LD 	Tip, 45-degree/Tip, straight: OM05510116/OM0551016 WITH 20g sleeve OPOS20L 	Tip, 45-degree/Tip, straight: OM05510116/OM0551016 WITH 19g sleeve OPOS19L 