

## SOVEREIGN Compact System With *WHITESTAR* Technology

**SOVEREIGN** Compact System provides access to **WHITESTAR** ICE Technology in a simple, user-friendly platform. The **SOVEREIGN** Compact System provides users with advanced fluidic performance, digital control over vacuum and power and patented Occlusion Mode for safer, more predictable surgery. With **SOVEREIGN** Compact System, lens removal procedures are performed in a closed system to prevent cross-contamination. **SOVEREIGN** Compact System has variable control over ramp speed and other performance parameters critical to maintaining a quiet, stable chamber.

The energy management capabilities of variable **WHITESTAR** ICE Technology provide surgeons with the ability to manage cutting power and efficiency to meet their individual preferences. **WHITESTAR** ICE Technology, with user friendliness and advanced fluidics, makes the **SOVEREIGN** Compact System an excellent option.

### NEW

Now available  
with advanced  
Ellips FX  
Technology!



**Phaco Handpiece**  
Part # SOV680290



**WHITESTAR**  
Phaco Handpiece  
Part # 690697



**ELLIPS FX**  
Handpiece  
Part # 690880



**Open Toe Foot Pedal**  
Part # CMP680720



**Closed Toe Foot Pedal**  
Part # CMP680700



**Remote Control**  
Part # CMP680135



**Single-Use Tubing**  
Part # OP061

**SOVEREIGN Compact**  
System Console  
Part # CMP680300

**SOVEREIGN Compact**  
System Cart  
Part # CMP680002

**SOVEREIGN Compact**  
DC Powered IV Pole  
Part # CMP680P2

Please ask your local AMO representative for the corresponding SOVEREIGN Compact owner's manual.