

tote-l-vac®

DESCRIPTION AND APPLICATION

Tote-L-Vac is a battery powered high vacuum and high flow suction unit. It features a soft, yet durable carrying case with easily accessible controls and battery. Although this unit is very compact, it provides high vacuum and high flow rates necessary for emergency oral aspiration procedures. Tote-L-Vac accepts both 800cc and 1200cc standard collection canisters. Tote-L-Vac is the perfect battery powered portable suction unit for EMS, home health, hospital patient transport and where ever battery power is a requirement.

FEATURES AND BENEFITS

- Vacuum regulator and gauge
- Housed in a durable carrying case
- Easy to clean and process for quick turnarounds
- Compact and lightweight

IDEAL FOR

It is an excellent choice for EMS, hospital patient transport, extended care facilities, doctors offices, home healthcare or where ever a battery powered suction unit is required.

BATTERY POWERED ASPIRATOR

E.M.S. • Patient Transport • Home Health Care • Surgicenters •
Nursing Homes • Physician's offices



TECHNICAL SPECIFICATIONS**Performance**

Vacuum Range: up to 550mmHg
 Free Air FLOW: 36 LPM
 Pump: 12 VDC oil-less diaphragm type

Collection Device

Canister: disposable plastic with mechanical shutoff
 Capacity: 800cc standard
 Tubing: 6' (182.7 cm) disposable sterile plastic
 Suction Tip: disposable Yankauer with control vent

Controls

Vacuum Regulator: rotary type on panel
 Vacuum Gauge: Calibrated to 0-30"/760 mmHg

Electrical Requirements

DC: 12 volts (nominal)

Battery:

Type: rechargeable sealed lead-acid
 Capacity: 12 volts, 3.0 ampere-hour
 Average Run Time: 30 minutes when fully charged
 Charger: fast/float charger
 UL/1310 listed
 Charge Time: 8 hours or less to 85% charging, depending on initial charge of battery
 Charge Method: 110 or 220 volt outlet, or 12 VDC vehicle outlet with optional cord

Physical

Weight: 7.0 Kg (15.5 lbs)
 Dimensions: 20.3 x 30.5 x 24.1 cm
 8"H x 12"W x 9.5"D

Warranty: One year

**PURCHASING SPECIFICATIONS**

tote-l-vac[®] Part Number: AI4900000 - 120 Volt
 749200 - 220 Volt

SHIPPING INFORMATION**Shipping weight:**

7.1 Kg (20 lbs.)

Package Size:

30.5 x 20.3 x 20.3 cm
 (12" H x 8" W x 8" D)

Packaging:

One per carton

AVAILABLE ACCESSORIES

| PART# | DESCRIPTION |
|-----------|---|
| AI4900400 | AC to DC Converter |
| AI5100125 | 12VDC Cigarette Power Cord |
| 749170 | Battery Charging Adaptor (for charging battery out of Tote-L-Vac) |
| 749146 | Replacement Lead Acid Battery |
| 749121 | AC Charger |
| AI4080010 | Disp. Collection Canister, 800cc (case of 10) |
| AI4080025 | Disp. Collection Canister, 800cc (case of 25) |
| 750150 | Filter,Hydrophobic with 5" Tubing (case of 50) |

